


## Tissue phantom for assessment of surface temperature of an ultrasonic transducer

### Schematic diagram of a typical setup to measure surface temperature



1. Ultrasonic transducer under test
2. Silicone skin-mimic (0.5 to 1.5 mm)
3. Agar tissue-mimicking material (TMM)
4. Thermal sensor (thermocouple)
5. Acoustic absorber

**Note:** Thermal sensor and acoustic absorber are not included as part of this phantom

### Specifications

The TMM complies with the specification given in clause 11 of IEC 60601-2-37 Ed 2, and is prepared following guidelines given in Annex DD of the same standard.

### Maintenance

The TMM should be stored in the closed container under normal laboratory conditions (18–25 °C). While it is stored, keep the TMM in a water/glycerol solution to prevent it from drying out and to avoid air contact. This solution contains 88.1 % (by weight) deionised water and 11.9 % (by weight) glycerol (purity > 99 %). While the TMM is stored, its surface should be submerged by at least 5 mm of solution.

The shelf life of the TMM, if it is preserved without air contact, is at least 1 year. To extend the shelf life to at least 2 years, add to the storage solution an anti-fungal agent: 0.5 % (by weight) benzalkonium chloride.

### Disposal

For disposal, local regulations should be followed. Alternatively, the samples may be returned (by prior arrangement) to the Acoustics group at NPL.

### Acoustic and thermal properties of materials (from Annex DD of IEC 60601-2-37 Ed 2)

Sample	Velocity C (m s <sup>-1</sup> )	Density ρ (kg m <sup>-3</sup> )	Attenuation coefficient α (dB cm <sup>-1</sup> MHz <sup>-1</sup> )	Acoustic Impedance Z (kg m <sup>-2</sup> s <sup>-1</sup> )	Spec. heat capacity C (J kg <sup>-1</sup> K <sup>-1</sup> )	Thermal conductivity κ (W m <sup>-1</sup> K <sup>-1</sup> )	Thermal diffusivity D (m <sup>2</sup> s <sup>-1</sup> )	Source of data
TMM	1540	1050	0.5	1.6 x 10 <sup>6</sup>	3800	0.58	0.15 x10 <sup>-6</sup>	TNO
Silicone rubber	1021	1243	1.8	1.3 x 10 <sup>6</sup>	-	0.25	-	TNO/Dow Corning

### Proportion of pure components in TMM

Component	% (by weight)
Glycerol	11.21
Water	82.95
Benzalkonium chloride	0.47
Silicon carbide (SiC (-400 mesh))	0.53
Aluminium oxide (Al <sub>2</sub> O <sub>3</sub> (0.3µm))	0.88
Aluminium oxide (Al <sub>2</sub> O <sub>3</sub> (3µm))	0.94
Agar	3.02
Sum	100