



深圳为尔康科技有限公司
 曾生：13632925349 QQ：274798107
 电话：0755-28896837
 网址：www.medicalqc.com
 地址：深圳市龙岗区沙平北路111号吉茂大厦608A

**UNMATCHED
 SOLUTIONS
 FOR QUALITY
 CONTROL**

Pro-Slit



**UPGRADED
 DESIGN**

The Pro-Slit phantom is the so called **slit camera** for **accurate measurement** of the **focal spot size according to IEC 60336:2005**. Its main advantages are **repeatability, accuracy and possibility to measure the size of any focal spot**. Because of the small attenuation of the outer cover this tool can also be used with dental digital detectors - detector receives enough radiation to be triggered.

**Technical data
 (can be modified to customer specifications):**

- dimensions: 40 x 40 x 11.5 mm
- 0.01 mm wide and 10 mm long slit with 8° spread
- diaphragm made of tungsten
- embedded in PMMA

Product features:

- complies with:
 - IEC 60336:2005
 - IEC 61223-3-1, 2, 4
- CE certified
- the Manual provides detailed guidelines for carrying out each test, results assessment and registration



深圳为尔康科技有限公司
 曾生：13632925349 QQ：274798107 电话：0755-28896837
 地址：深圳市龙岗区沙平北路111号吉茂大厦608A 网址：www.medicalqc.com