

OxSim Flex® Specifications

Physical Dimensions			
Size	2.5 in x 2.5 in x 1.25 in (6.4 cm x 6.4 cm x 3.8 cm) excluding finger		
Weight	5 oz (142 g) (excluding battery)		
SpO2 Connection	Optical		
Environmental			
Operating	Temperature: 10°C to 35°C (50°F to 95°F), Relative Humidity: 10 to 80% non-condensing, Altitude: 3,000m (9,843ft)		
Storage	-20° C to $+60^{\circ}$ C (-4° F to $+140^{\circ}$ F)		
Power			
Battery Power	1 AA battery (any chemistry)		
Battery Life	 Continuous: OxSim Flex: 10+ hours, OxSim Flex with Bluetooth: 6.5 hours Typical Use (Non Continuous): OxSim Flex: 20+ hours, OxSim Flex with Bluetooth: 13 hours 60 days in sleep mode (unit sleeps after 10 min. with no probe connection) 		
AC Operation	Included Class II limited power source 6VDC supply operates from 100-240VAC, 50-60Hz		
Saturation	•		
Range	10% – 100%		
Accuracy	Saturation 85% - 100% 70% - 84% 60% - 69%	Autoprobe Selected ± 2% of reading ± 3% of reading NA	Manufacturer Probe Selected ± 1% of reading ± 2% of reading ± 3% of reading
Heart Rate			
Range	25 – 240 BPM		
Accuracy	± 2 BPM		
Perfusion Index	-		
Range	10% - 100% (0.2 – 2.0 PI)		



OX-2 OxSim Flex® Optical SpO2 Pulse Oximeter Simulator Features:

- Complete control of simulations!
- Saturations in 1% increments (10-100%) •
- Pulse rates in 1 BPM increments (25-240)
- Create custom preset simulations for faster Alarm Testing •
- Perfusion Index levels in 1% increments (10-100%) •
- Factory Preset Simulations: •
 - Saturation: 85%, Heart Rate: 80 bpm, Perfusion Index: approx. 2.0 0
 - Saturation: 95%, Heart Rate: 40 bpm, Perfusion Index: approx. 2.0 0
 - Saturation: 98%, Heart Rate: 80 bpm, Perfusion Index: approx. 2.0 0
 - Saturation: 98%, Heart Rate: 140 bpm, Perfusion Index: approx. 2.0 0
 - Saturation: 99%, Heart Rate: 80 bpm, Perfusion Index: approx. 2.0 0
- Masimo Rainbow SET® Compatible •
- Automatic detection of SpO2 Manufacturer •
- Color Display OLED Technology •
- Only 5 oz. (0.14 kg, excluding battery) •
- Fast and Accurate •
- Battery Operated minimum 10 hours •
- Power Save Mode •
- Rugged: Polycarbonate housing. Drop tested 65 times from a height of 3 feet (1m) to a hard • surface with no failures!
- Field-upgradable Software •
- Four-year Limited Warranty

Package Contents:

- (1) OxSim[®] Flex[®] Optical SpO2 Pulse Oximeter Tester •
- (1) Power Supply (6VDC) •
- (1) AA battery •
- (1) Certificate of Calibration •
- (1) Carrying Case



深圳为尔康科技有限公司 曾生:13632925349 QQ:274798107 MedicalQC 电话: 0755-28896837 地址:深圳市龙岗区沙平北路111号吉茂大厦608A 网址:www.medicalgc.com