

UMR-3D



The UMR-3D is the power meter you need to service the 1 and 3 MHz continuous output therapy level ultrasound generators.

UMR-3D power meter offers a big value in a portable package designed to match your service requirements and priced to fit your budget.

Features:

- Accepts 1 and 3 MHz transducers
- Transducer holder
- Quick response LCD
- Uses 1 oz. degassed water for coupling
- Frequency selector switch
- Fully contained in a compact case

SPECIFICATIONS

Input Power:	0-30 Watts	Operating Temperature:	60°F - 80°F
Resolution:	0.1 Watts	Weight:	8 lb. / 3.8 kg
Accuracy:	±3%	Coupling:	1 oz. Degassed Water
Linearity:	0.5%	Power Supply:	120V, 60 Hz
Options:	220V 50 Hz operation RS-232C communication port with Test Software	Case:	6" x 9" x 6" / 150 x 230 x 150 mm
		Applicable Transducer:	10 cm² or 5 cm²

Notes:

For your convenience 1 oz. Degassed Coupling Water can be provided in packages of twelve. Annual recalibration is recommended.



深圳为尔康科技有限公司 联系人：曾祥满 手机：13632925349

QQ：274798107 电话：0755-28896837 地址：深圳市龙岗区沙平北路111号608A

网址：www.medicalQC.com 邮箱：szchina1718@163.com